

High-class portable lung ventilator with a built-in flow generator and modern touch screen management.

The screen displays all digital parameters, as well as respiratory graphics, including FiO2, EtCO2 monitoring in real-time mode.

The ventilator is applicable for adults and pediatric patients

weighting at least 5 kg.

JSC "Ural instrument-engineering plant"

the 25th km of the Chelyabinsk path, Sysertsky area, Sverdlovsk region, Russia 620 000 tel: +7 343 359 93 15; fax: +7 343 359 93 15; e-mail: e.zolotnitskaya@upz.ru; www.upz.ru





New High-class lung ventilator

"AVENTA-M" New opportunities for

autoadaptive control





High-class lung ventilator "AVENTA-M"

High-class portable lung ventilator with a built-in flow generator (turbine) and modern touch screen management. The screen displays all digital parameters, as well as respiratory graphics. The ventilator is applicable for newborns 3kg, children and adults.

Features:

- High sensitive correction of the ventilation in case of independent breathing efforts of the patient;
- Automatic adjustment of ventilation parameters and alerts limits with due account for patient's weight;
- Maneuver of therapeutic bronchoscopy;
- Automatic selection of pressure increasing optimal rate (Pramp AUTO);
- Sensitive smart trigger does not miss real breathing attempts and minimizes the probability of false
- alarm even under conditions of leaks and NIV ventilation;
 Automatically controlled Exhalation trigger (Esens-Auto);
- Automatically controlled FiO2;
- Automatically controlled PEEP;
- Auto maintenance (support) of set Minute Ventilation in all modes;
- Automatic leakage compensation;
- Expiratory valve (unit) installed inside the lung ventilator;
- 5 hours internal battery;
- Ventilation pattern adjustment to the individual lungs characteristics (AUTO-MVG);
- EtCO2 and FiO2 monitoring;
- Wireless data transfer;
- Ventilation map (display of relation between Minute ventilation and breath rate, FiO2 level and PEEP rate, Spontaneous breath rate (f_spont) and Inspiratory pressure (Pi);
- Multiscreen (simultaneous display of real lung's image, Paw, Flow and VT graphics);
- Rapid shallow breath index (RSBI);
- Standards rapid shallow breath index (RSBN). Standards RSBN index of the patient's ideal body weight (IBW)

Perfomance specifications:

Breath rate (f)	1-150 bpm
SIMV breath rate	
Respiratory volume (VT)	10-3000 ml
Inhalation Peak flow (Peak flow)	180 l/min
Plateau time(Tpl)	0-2c
PEEP (including newborns ≥5kg)	0-50 с м Н ₂ О
CPAP	
FiO2	21 – 100%
Inspiratory pressure (Pi)	0-100cмH ₂ O
Threshold pressure drop	
Inspiratory time (Ti)	0.2 -15 c
I : E Ratio	(1:10) - (4:1)
Pressure support(PS)	0-100 смH ₂ O
Flow-trigger (F-trig)	0.5-25 l/min, OFF
Pressure trigger (P-trig)	0.5-25cмH ₂ O, OFF
Inspiratory pause	
Expiratory pause	
FiO2 increasing	(100% O ₂ for 2 min)
Color TFT touch-screen display 12 inches	
Block of touch panel	
WOB_tot Shows the total work performed by the	e device and the patient on inhalation

WOB_pat Shows the work performed by the patient on inhalation WOB_imp Shows the work performed by the patient on inhalation triggering (overcoming the breathing circuit resistance)	
_ , , , ,	
(0.0.0000000)	
Pause time 15-60 sec	
Time of pressure increasing (Pramp) 50-500mc, AUTO	
Built-in nebulizer providing < 3 microns and can be used in both TURNED OFF and TURNED ON positions	
SpO2 and PR	

Monitoring parameters:

 Peak Inspiratory pressure (Ppeak) Exhalated Tidal Volume (Vte) Mean Airway Pressure (Pmean) Inspirated Tidal Volume (Vti) Minute volume (MV, Mv_spont, MV_mand)
 Plateau pressure (Pplat) Positive End Expiratory Pressure (PEEP) • Breath rate (f, f_spont, f_mand) • FiO2 Inspiratory time (Ti) Compliance (Cstat, Cdyn) Expiratory time (Te) Resistance (R stat) • I :E Ratio Summery frequency Spontaneous minute control Trends (changed parameters) Inhalation time constant (Tc i) of up to 50 parameters Exhalation time constant (Tc e) with 60 seconds resolution End exhalation pressure (Pe end)

Parameters displayed

Accumulator capacity, patient data, alarms, settings, clock, modes (regimes) etc.

Graphic system displays the following graphics:

- · Waves (Stream, Pressure, Volume);
- loops (Flow/Volume, Pressure/Volume, Flow/Pressure) trends (up to 7 days)

Gas source

- compressed gas distribution system, oxygen Concentrator, balloon (O₂)
- internal turbine (AIR)

Ventilation modes:

A/C-VC (controlled volume) A/C-PC (controlled pressure) SIMV-VC SIMV-PC SIMV-PS (pressure support) SIMV-P/S-V/G CMV **PEEP** CPAP (Spont.), CPAP/PS CPAP/PS-VG Dual-Level (biphase breathing) Dual-Level/PS-VG NIV (noninvasive ventilation) Auto MVG Apnoe /additional ventilation **APRV**